

Details for AQB Implant tools, sets and separately-sold parts

●Settools

Product name	Model	Manufacturer	Authorization/certification/ registration number
AQB Round bur	-	Dentech Corporation	21900BZX00404000
AQB guide drill	S L		
AQB direction checker	S	Dentech Corporation	13B2X00094000072
	M L		
AQB spiral drill	φ2.5S	Dentech Corporation	21900BZX00404000
	φ2.5L		
	φ3.5S		
	φ3.5L		
	φ4.5S		
	φ4.5L		
AQB counter drill	φ3.5 φ4.5		21900BZX00404000
AQB end mill reamer	φ2.75 φ3.75		13B2X00094000065
	φ4.75		
AQB reamer extension	-		13B2X00094000075
AQB finger driver	S		13B2X00094000073
	L		
AQB free wrench	-		13B2X00094000070
AQB free wrench adapter	-		13B2X00094000074
AQB alarm gauge	For φ3	AQB-ABI Implant Co.,Ltd. Dentech Corporation	11B1X10027080901 13B2X00094000071
	For φ4		
	For φ5		
AQB reamer	φ2.9	Dentech Corporation	13B2X00094000066
	φ3.9		
	φ4.9		
AQB circular knife	For φ3		21900BZX00404000
	For φ4		
	For φ5		
AQB fixer	For φ3		13B2X00094000067
	For φ4		
	For φ5		
AQB fixersupport	-		13B2X00094000068
AQB fixerjoint	-		13B2X00094000069
AQB fixerextension	-		13B2X00094000076
Tool case	-	-	-
X-ray guide	-	-	-
Tool guide (set of 6 for φ3, φ4 and φ5)	-	-	-

●One-piece T-type dedicated tool

Product name	Model	Manufacturer	Authorization/certification/ registration number
AQB fixer T-type	For φ3	Dentech Corporation	13B2X00094000067
	For φ4		
	For φ5		
AQB counterbore drill	For φ3	Dentech Corporation	21900BZX00404000
	For φ4		
	For φ5		

AQB・ABI Implant Co., Ltd.
Head Office: 899-5 Kitahara-cho, Tokorozawa City, Saitama Prefecture, Japan 359-0004

Orders & inquiries Higashi-Nihonbashi 1-chome Building 8F, 1-4-6 Higashi-Nihonbashi, Chuo-ku, Tokyo, Japan 103-0004
Dental Business Tel: +81-3-5839-2541 Fax: +81-3-3862-1264 Email: info@aqb.jp Website: https://www.aqb.jp
Headquarters

*This list,as of November 2024,is subject to change without notice.

●Two-piece dedicated parts

Product name	Model	Manufacturer	Authorization/certification/ registration number
Abutment screw	-	Advance Co., Ltd. AQB-ABI Implant Co.,Ltd.	20500BZZ01172000

●Two-piece dedicated tools

Product name	Model	Manufacturer	Authorization/certification/ registration number
HEX driver bit *1)	φ3.0S	Dentech Corporation	21900BZX00404000
	φ2.35S		
	φ2.35L		
AQB bone profiler	4		
	5		
AQB bone profiler guide	-		
SOL driver bit	φ3.0S	AQB-ABI Implant Co.,Ltd.	11B1X10027060404
Torque control	-	Yoshida Dental Mfg. Co., Ltd.	13B1X00133Z00028
HEX driver mini	-	Dentech Corporation	13B2X00094000126
AQB custom type	4	AQB-ABI Implant Co.,Ltd.	11B1X10027000031
	5		
Abutment replica *2) (21 types x 3 pieces)	-	Advance Co., Ltd. AQB-ABI Implant Co.,Ltd.	13B1X00205000030 11B1X10027000030
*1) For finger driver: φ3.0S. For torque control: φ2.35S, φ2.35L *2) Replica single item (5 pieces of each size)			

●Two-piece T-type dedicated tools

Product name	Model	Manufacturer	Authorization/certification/ registration number
AQB counterbore drill	(2P) φ4-5	Dentech Corporation	21900BZX00404000

●Optional tools

Product name	Model	Manufacturer	Authorization/certification/ registration number
AQB reaming drill	φ2.75	Dentech Corporation	21900BZX00404000
	φ3.75		
	φ4.75		
	φ2.75L		
	φ3.75L		
	φ4.75L		
Surgical bur	164RF/018	Sato Dental Co., Ltd.	13B2X00094000120
AQB torque wrench	-		13B1X10122000009

●Micro motors

Product name	Model	Manufacturer	Authorization/certification/ registration number
Surgic Pro2 (with light)	LED	Nakanishi Co., Ltd.	303ALBZX00014000
Surgic Pro2 (without light)	No LED		

●Other

Product name	Model	Manufacturer	Authorization/certification/ registration number
Leaflet (100 copies)	-	AQB-ABI Implant Co.,Ltd.	-
Treatment explanation leaflet (100 copies)	-		-



One-piece type												Brand: AQB Implant Auth. number: 20500BZZ01172000											
φ 3	Model No.	3SS		3SM		3SL		3MS		3MM		3ML		3LS		3LM		3LL		3Y9		3Y11	
	11																						
	9																						
	7																						
mm	0																						
	8																						
	10																						
	12																						
mm	14																						
φ 4	Model No.	466		468		4SS		4SM		4SL		4MS		4MM		4ML		4LS		4LM		4LL	
	11																						
	9																						
	7																						
mm	0																						
	6																						
	8																						
	10																						
mm	12																						
φ 5	Model No.	566		568		5SS		5SM		5SL		5MS		5MM		5ML		5LS		5LM		5LL	
	11																						
	9																						
	7																						
mm	0																						
	6																						
	8																						
	10																						
mm	12																						

: Delivery may take 1-2 months after the order due to the special size.

One-piece T-type												Brand: AQB Implant Auth. number: 20500BZZ01172000				
φ 3	Model No.	T3SS		T3SM	T3SL		T3MS	T3MM	T3ML	T3LS		T3LM	T3LL	T3Y9		T3Y11
	11															
	9															
	7															
mm	0															
	8															
	10															
	12															
mm	14															
φ 4	Model No.	T46S		T468	T4SS		T4SM	T4SL		T4MS	T4MM	T4ML	T4LS		T4LM	T4LL
	11															
	9															
	7															
mm	0															
	6															
	8															
	10															
mm	12															
φ 5	Model No.	T56S		T568	T5SS		T5SM	T5SL		T5MS	T5MM	T5ML	T5LS		T5LM	T5LL
	11															
	9															
	7															
mm	0															
	6															
	8															
	10															
mm	12															

: Delivery may take 1-2 months after the order due to the special size.

Indication		One-piece type			One-piece T-type		
	Outer diameter	4			T4		
		M			M		
	Coated part	6 : 6mm 8 : 8mm 10 : 10mm 12 : 12mm 14 : 14mm			6 : 6mm 8 : 8mm 10 : 10mm 12 : 12mm 14 : 14mm		
		Abutment			Abutment		
	Abutment	3 : 3mm 4 : 4mm 5 : 5mm			T3 : 3.4mm T4 : 4.4mm T5 : 5.4mm		
		6 : 6mm 7 : 7mm 8 : 8mm 9 : 9mm 10 : 10mm 11 : 11mm			S : 7mm 8 : 8mm 9 : 9mm 10 : 10mm 11 : 11mm		

AQB Implant two-piece system combination

The diagram illustrates the AQB Implant two-piece system combination, showing various components and their combinations. The components are organized into several categories:

- Healing cap**: 4, 5
- Healing abutment**: SS, S, M, L
- Coping screw**
- AQB pick-up coping**: $\phi 4$, $\phi 5$
- Multi abutment**: $\phi 4$, $\phi 5$
- Abutment screw**
- Straight abutment**: 5-6-0, SS, S, M, L
- Straight abutment W**: WS, WM, WL
- Angle abutment**: 15°SS, 15°S, 15°M, 15°L, 25°SS, 25°S, 25°M, 25°L
- Angle abutment W**: 15°WS, 15°WM, 15°WL, 25°WS, 25°WM, 25°WL
- Multi abutment**: $\phi 4$, $\phi 5$
- AQB custom type**: 4, 5
- OD abutment**: A, B
- OD screw**
- OD abutment base**
- Implant analog**: $\phi 4$, $\phi 5$
- Two-piece fixture**: 4081, 4101, 4121, 4082, 4102, 4122, 5081, 5101, 5121, 5082, 5102, 5122
- Two-piece T-type fixture**: T4121, T4141, T4122, T4142, T4123, T4143, T4125, T4145


The diagram also includes a vertical arrow indicating the flow of components from the top (Healing cap) to the bottom (Implant analog).

Common use of abutments			
		φ4 fixture	φ5 fixture
		4A or 4B (left figure)	5A or 5B (left figure)
Straight abutment Angle abutment	5-6-0	0.5mm	0.0mm
	SS	0.5mm	0.0mm
	S	1.0mm	0.5mm
	M	2.0mm	1.5mm
	L	3.0mm	2.5mm
	WS	1.5mm	1.0mm
Healing abutment	WM	2.5mm	2.0mm
	WL	3.5mm	3.0mm
	SS	1.0mm	0.5mm
	S	3.0mm	2.5mm
Multi abutment	M	4.0mm	3.5mm
	L	5.0mm	4.5mm
	4	4.0mm	-
OD abutment base	5	4.5mm	4.0mm
	-	1.5mm	1.0mm

*Abutment screws are also sold separately.

Healing abutment

Brand: AQB Implant
Auth. number: 20500BZZ01172000



Indication


Gingiva part length
5mm

M

Gingiva part length
SS : 0.5mm
S : 2.5mm
M : 3.5mm
L : 4.5mm

Multi abutment

Brand: AQB Implant
Auth. number: 20500BZZ01172000




Indication

Outer diameter
4 : 4mm
5 : 5mm

OD abutment base

Brand: AQB Implant
Auth. number: 20500BZZ01172000




Indication

1mm
5mm

AQB pick up coping

Brand: AQB pick up coping
Auth. number: 11B1X10027060402




Used with a multi abutment when taking an impression.

Indication

Inner diameter
4 : 4mm
5 : 5mm
Outer diameter
4 : 7mm
5 : 8mm

AQB Implant analog

Brand: AQB Implant analog
Auth. number: 11B1X10027060401




Indication

Outer diameter
4 : 4mm
5 : 5mm

AQB coping screw

Brand: AQB coping screw
Auth. number: 11B1X10027060405




Used with a multi abutment and AQB pick-up coping when taking an open tray impression.

Indication

Outer diameter
2.5mm
11.5mm
17.8mm

OD abutment (burn-out type)

Brand: OD Abutment
Auth. number: 13B1X00205060403
11B1X10027060403



For use with the OD Abutment base.

Indication



Model No.
A
B

Outer diameter
5.04mm
Inner diameter
4.1mm
2.6mm
4.2mm
5mm
1.8mm

OD Abutment A OD Abutment B

2-piece fixture

Brand: AQB Implant
Auth. number: 20500BZZ01172000




Collar length
Full length
Outer diameter
Outer diameter 2

Indication
4 1 0 2
Outer diameter 1
4 : 4mm
5 : 5mm
Full length
0.8 : 0.8mm
10 : 10mm
12 : 12mm
Collar length
1 : 1mm
2 : 2mm

*The collar is not coated

Two-piece fixture T-type

Brand: AQB Implant
Auth. number: 30100BZX00252000



Collar length
Outer diameter
Full length
Outer diameter 1
Outer diameter 2


Indication
T4 1 2 5
Outer diameter 1
T4 : 5mm
Outer diameter 2
T4 : 4mm
Full length
12 : 12mm
14 : 14mm
Collar length
1 : 1mm
2 : 2mm
3 : 3mm
5 : 5mm

*The collar is not coated

*Abutment screws are also sold separately.

Straight abutment

Brand: AQB Implant
Auth. number: 20500BZZ01172000



Indication


Gingiva section
3mm
0mm
6mm
5mm

Gingiva section
3mm
0mm
6mm
5mm

Indication
M
5-6-0 : 0 mm
SS : 0 mm
S : 0.5mm
M : 1.5mm
L : 2.5mm

Straight abutment W

Brand: AQB Implant
Auth. number: 20500BZZ01172000




Indication

Gingiva section
6mm
5.9mm

Indication
WM
WS : 1mm
WM : 2mm
WL : 3mm

Angle abutment

Brand: AQB Implant
Auth. number: 20500BZZ01172000



Indication


Gingiva section
6mm
5mm

Angle
25°
15°
25°

Indication
M
SS : 0 mm
S : 0.5mm
M : 1.5mm
L : 2.5mm

Angle abutment W

Brand: AQB Implant
Auth. number: 20500BZZ01172000



Indication

Gingiva section
6mm
5.9mm


Angle
25°
15°
25°

Indication
WM
WS : 1mm
WM : 2mm
WL : 3mm

*1: The 15° (SS, S, M, L) abutment has two types: one for use with 4 fixtures only, and one for use with 4 and 5 fixtures. The 25° (SS, S, M, L) abutment has one type: for use with 4 and 5 fixtures. The 4 abutment is a special size so that the delivery may take 1-2 months after order.

Procedure

Tools used (common)




Related products


Tool set



Tool guide



X-ray guide



Abutment replica set



Surgic Pro2

